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
Product Specifications

 DS-TDSB00-EKT/4m	Working Frequency	60 to 64 GHz	Horizontal & Vertical FoV	-45° to 45°, -45° to 45°
	Modulation Wave	FMCW	Communication Mode	Wi-Fi, RS-485 interface
	Frequency Span	2 GHz	Working Voltage	9 to 12 VDC
	Detection Range	0.1 to 4 m	Working Electric Current	≤ 200 mA @ 12 VDC
	Data Cycle	70 ms	Weight	133 g

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Product Specifications

 DS-TDSB00-EKH/2m	Working Frequency	60 to 64 GHz	Horizontal & Vertical FoV	-45° to 45°, -45° to 45°
	Modulation Wave	FMCW	Breath/Heart Rate Resolution	0.08 Hz
	Frequency Span	3.5 GHz	Data Cycle	200 ms
	Detection Range	0.1 to 2.7 m	Communication Mode	Wi-Fi, RS-485 interface
	Weight	125 g	Working Voltage & Electric Current	9 to 12 VDC, ≤ 200 mA @ 12 VDC

SOLUTION

Hikvision Radar-based Elder Care Solution- Protecting the safety of the elderly and improving the operation efficiency of nursing homes

——Charlie.Zhou, zhouzihao@hikvision.com

Introduction to the Elderly Care Industry

1. The era of global aging has begun

According to official statistics, by the end of 2021, the global population over the age of 65 has exceeded 700 million, accounting for about 9% of the total population. Among them, the aging situation in developed countries and regions is more serious. In Japan, for example, the proportion of aging has reached 29%, and in Europe, it has reached 20.3%. It means that about one in four people is the elderly, and aging population will continue to increase, the elderly population is likely to double by 2050.

2. Nursing home beds and nursing staff are in short supply

The direct impact of aging is that the construction speed of nursing homes cannot meet the rapid growth of the elderly population. The preparation and construction of a nursing home takes a long time, especially in institutions with better environmental conditions, beds are seriously scarce, and the elderly even need to queue up until 2-3 years later to move in.

At the same time, the shortage of nursing staff in nursing homes is also a big issue. Because caregivers are in a relatively poor situation, they usually have to work long hours, but with low incomes. In addition, caregivers often need very professional skills and there is often a high turnover rate. In the UK, job vacancies for long-term care workers reached 105,000, and the 8.2% vacancy rate was even more than double of the average job vacancy in the UK.

3. Digital development of nursing homes

At present, the traditional nursing home has been unable to meet the market demand, in the case of ensuring the basic life safety of the elderly, the nursing home needs to pay more attention to the spiritual needs of the elderly. The construction of nursing home has also developed, such as libraries, movie theaters, universities for the elderly, activity centers also become the next target of nursing homes.

Therefore, the digital upgrading of nursing homes is particularly important, and the access of intelligent equipment can greatly improve the management and operation efficiency of the organization. Devices such as fall monitoring and heart rate monitoring can provide more health supervision for the elderly. Personnel positioning and behavior analysis can give the elderly more space for activities and freedom.

Solution Overview

Around the elderly health management, we offer two radar products to guard and detect the health of the elderly in the indoor scenarios: Fall Detection Radar and Auxiliary Care Radar.

1. Fall Detection Radar

With the increase of age, the physiological functions of the elderly are seriously degraded, and it is easy to fall or have other accidents. Falls are the leading cause of injury death among the elderly over the age of 65. Even in nursing homes, about 30 percent of the elderly fall every year. Among them, more than half of serious accidents are caused by not being detected and dealt with in time after a fall. The Fall Detection Radar is used to detect the fall accident of the elderly in the indoor scene, and send an alarm to the nurse station to shorten the rescue time.

1) No privacy disclosure

In the interior scene, bedrooms and bathrooms are the main areas for elderly people's activities. Among them, uneven floors, slippery floors, and wet bathrooms are prone to fall down. In addition, sharp bed corners and table corners are also difficult to avoid bumps. These falls are very fatal for the elderly. For such private spaces, traditional video monitoring cannot meet such scenarios. Hikvision Fall Detection Radar solution uses millimeter wave radar cloud dot matrix imaging technology and captures the height difference caused by the fall of the elderly, combined with posture recognition, can effectively protect user privacy, provide high accuracy of fall judgment, and bring more humanized services to users.

2) Non-contact


In the event of a serious fall, the elderly are usually unable to call or ask for help. The radar can actively detect and identify personnel changes, and will send an alarm to the nurse center after the fall, which greatly shortens the time from accident to discovery and avoids secondary injuries caused by delayed rescue.

3) Integration with third-party platforms

HikCentral Professional provides holistic management including device management, data analysis, and real-time alarm notifications. At the same time, data can be integrated with 3rd party platforms through OpenAPI.

2. Auxiliary Care Radar

Nursing homes, as the place where the elderly live for a long time, need to provide 24-hour caring and supervision, and even during the night, the emergencies of the elderly are uninterrupted. Manual inspection is labor-intensive and inefficient. **Based on such problems, intelligent auxiliary nursing equipment can reduce costs and provide effective safety supervision.**

1) In bed signs & body movement detection

When the elderly are in bed, the radar can detect the heart rate, breathing rate, and record physical movement data such as rolling over and getting up, which will be clearly displayed on the platform. Long-term physical vital signs data can be used to assist in judging the health status of the elderly.

2) Sleep detection

The Auxiliary Care Radar can analyze and record the sleep quality of the elderly (deep sleep, light sleep, and awake) through its own algorithm, and the nursing home can view abnormal sleep reports to achieve refined nursing.

3) Abnormal leave bed alarm

Real-time monitoring of the status of the elderly, including whether to return to the room at night, whether to get up on time in the morning, etc. Caregivers do not have to frequently inspect the room to check the status of the elderly which will greatly improve management efficiency. At the same time, the radar can provide an abnormal status alarm, such as the elderly getting up at night or leaving the bed for too long

because of a come or a fall. Or maybe they didn't get up at the right time because they weren't feeling well. After detecting these abnormalities, the radar will send out a real-time alarm to remind the nurse to check the situation, providing assistance as soon as possible after an accident.

4) Single/multiple-person room management

In addition to the single room mode, the multi-room mode is also supported. Each bed corresponds to a radar, and the body and sleep indicators corresponding to the elderly in each bed can be clearly displayed in the platform, as well as the fall alarm and emergency rescue that occur in the room, so as to realize the refined management and operation of the nursing home. The platform enables the export of relevant signs and sleep data for use by third-party professional organizations.

5) Nursing home staff preview

The status of the elderly in bed is displayed in a list, which is convenient for nurses to identify issues at a early stage, greatly improving the efficiency of room inspection. At the same time, by looking at the real-time heart rate and breathing rate of the elderly in bed, caregivers can detect abnormalities in the physical condition of the elderly early.


The screenshot displays the HikCentral Professional Web Client interface. At the top, there are navigation options: Room Management, Export, and Update. A search bar for rooms is also present. Below the navigation, there are several status indicators: Leave Room (24 case), Leave Bed (6 case), SOS (1 Time), Fall (1 Time), Intruder Alarm (0 Time), and Abnormal Leave Bed Alarm (0 Time). The main area features a table with the following data:

Room	No	Name	In the room	In the bed	Current Heart Rate	Current Breath Rate
302	01	xxxxxx	Y	N	70	16
302	02	xxxxxx	N	Y	-	-
302	03	xxxxxx	Y	N	-	-

At the bottom, there is a pagination control showing 'Total 100' and a page number '1' out of '100'.

3. System topology

- 1) The radar supports access to the PC terminal and APP terminal of HikCentral Professional, and triggers an alarm after events such as falls and abnormal bed leaving.
- 2) For some nursing homes with their own operating platforms, they can also connect the radar to their own platforms through OpenAPI for equipment management and alarm integration.

Core Products

1. Fall Detection Radar

Product Specifications

 DS-TDSB00-EKT/4m	Working Frequency	60 to 64 GHz	Horizontal & Vertical FoV	-45° to 45°, -45° to 45°
	Modulation Wave	FMCW	Communication Mode	Wi-Fi, RS-485 interface
	Frequency Span	2 GHz	Working Voltage	9 to 12 VDC
	Detection Range	0.1 to 4 m	Working Electric Current	≤ 200 mA @ 12 VDC
	Data Cycle	70 ms	Weight	133 g

2. Auxiliary Care Radar

Product Specifications

 DS-TDSB00-EKH/2m	Working Frequency	60 to 64 GHz	Horizontal & Vertical FoV	-45° to 45°, -45° to 45°
	Modulation Wave	FMCW	Breath/Heart Rate Resolution	0.08 Hz
	Frequency Span	3.5 GHz	Data Cycle	200 ms
	Detection Range	0.1 to 2.7 m	Communication Mode	Wi-Fi, RS-485 interface
	Weight	125 g	Working Voltage & Electric Current	9 to 12 VDC, ≤ 200 mA @ 12 VDC